

# PEDAL OPERATED SUCTION PUMP

The Clements Pedal Operated suction pump is a robust portable, manually operated suction pump with vacuum indication and fluid collection, providing suction where electricity is not available. The Clements Pedal Operated Obstetric Suction System is also available to assist in obstetric deliveries and pharyngeal suction of the mother.

Vacuum storage feature removes need to continuously operate pump.

Metal handpiece is normally occluded and opened with thumb-operated spring-return release valve.

- Maximum Vacuum: -80 kPa [ -600 mmHg]
- Maximum Flow : 32 litres/minute
- Independent of battery and mains power
- Robust construction



# PEDAL OPERATED SUCTION PUMP

TECHNICAL SPECIFICATION		
	SUC 82010	SUC 82012
Standard Accessories	SUC 80350 - Metal Handpiece	SUC 80150 - Set of 3 Obstetric Cups
Maximum Vacuum	-80 kPa, [-600 mmHg]	
Maximum Flow Rate	32 L/min	
Gauge	63 mm diameter, Bourdon tube type, dual scale, graduated in kPa and mmHg	
Gauge Range	0 to -100 kPa graduated at 5 kPa 0 to -760 mmHg graduated at 20 mmHg	
Pump Type	Piston	
Manual Requirements	About 20 consecutive strokes to achieve maximum vacuum	
Collection Canister - Reusable	2 Litre autoclavable, shatter-resistant polycarbonate jar	
Overflow Protection - Reusable Canister	Float valve mechanism	
Filter	Disposable 3 micron anti-bacteria filter at inlet	
Dimensions	335H x 185W x 320D mm	
Weight	3.8 kg	
GMDN	47368	
Device Class	Class I	
Operation Mode	Continuous operation	

PART NUMBERS	
SUC 80150	Set of 3 Obstetric Cups, 40, 50, 60 mm with Puller Chain and Valve Assembly
SUC 80151	40 mm Obstetric Cup,with Puller Chain and Valve Assembly
SUC 80152	50 mm Obstetric Cup,with Puller Chain and Valve Assembly
SUC 80153	60 mm Obstetric Cup,with Puller Chain and Valve Assembly
SUC 80297 001	Yellow Medical Suction Tubing, 20 m
SUC 80301	2 L Reusable Collection Jar (Jar Only)
SUC 80312	2 L Reusable Collection Jar Complete with Bung Assembly and Overflow Cut-Off Valve
SUC 80350	Metal Handpiece, Finger Release Valve and Multi-Hole Tip
SUC 84200 033	Bung Assembly and Overflow Cut-Off Valve for 2L Collection Jar

All specifications are subject to change without notice. Please check the exact specification when ordering.

BMDi TUTA Healthcare Limited  
Unit 4B, 128-130 Frances St  
Lidcombe NSW 2141, Australia  
Telephone +61 2 9466 5300  
Email [info@medaust.com](mailto:info@medaust.com)  
Web [www.medaust.com](http://www.medaust.com)

Rev 5. © Copyright 2017 Medical Australia Limited.

**CLEMENTS**