

ACEEVAC PORTABLE SUCTION PUMP

MAINS OPERATED

ACeevac continues the tradition of the Clements Emergency pump, providing the same dependable level of performance in a lighter, more compact form factor.

Designed for high vacuum performance, reliability and minimal maintenance, ACeevac is ideal for clinical and home use or as an alternative to wall suction.

- Maximum Vacuum: -80 kPa [-600 mmHg]
- Maximum Flow : 21 litres/minute
- Direct vacuum regulation
- Compact and lightweight



ACEEVAC PORTABLE SUCTION PUMP

TECHNICAL SPECIFICATION	
Power Rating	220 - 240V 50/60Hz 110 VA 0.5A
Vacuum	-20 to -80 kPa, [-150 to -600 mmHg]
Maximum Flow Rate	21 L/min (measured through 1L reusable jar) (36 L/min measured at pump)
Gauge	40 mm diameter, Bourdon tube type, dual scale, graduated in kPa and mmHg
Gauge Range	0 to -100 kPa graduated at 5 kPa 0 to -760 mmHg graduated at 20 mmHg
Vacuum Regulator	Needle valve type
Pump Type	Piston - motor and pump combined
Motor Type	PSC AC motor
Collection Canister - Reusable	1 Litre autoclavable, shatterproof polycarbonate jar
Overflow Protection - Reusable Canister	Float valve mechanism / hydrophobic filter
Collection Canister - Disposable Liner	1 Litre disposable liner support polycarbonate jar (Optional)
Overflow Protection - Disposable Liner	Inbuilt hydrophobic and antibacterial filter (Optional)
Power Indicator	Neon indicator shows power available
Dimensions	360L x 215W x 195H mm
Weight	3.6 kg
GMDN	36777
Device Class	Class IIA
Electrical Protection	Class II
Operation Mode	Continuous operation

PART NUMBERS	
SUC 81025	High Vacuum High Flow Portable Suction Pump with 1 x 1L Reusable Collection Jar
SUC 81030 001	Kit, 1L Reusable Collection Jar (Pack 1) (1L Reusable Jar, Lid)
SUC 81030 002	Kit, 1L Disposable Liner Collection Jar (Pack 1) (1L Support Jar, 1L Disposable Liner)
SUC 81030 121	Kit, 1L Disposable Liner (Pack 50) (1L Disposable Liner)
SUC 81030 012	Kit, Hydrophobic and Antibacterial Inline Filter (Pack 10)

All specifications are subject to change without notice. Please check the exact specification when ordering.



BMDi TUTA Healthcare Limited
Unit 4B, 128-130 Frances St, Lidcombe
NSW 2141, Australia
Telephone +61 2 9466 5300
Email info@medaust.com

Rev 3. © Copyright 2017 Medical Australia Limited.

CLEMENTS